



AQUILA MEDICAL

天鹰医疗 产品介绍

Aquila Medical Product Introduction

天鹰刨削系列 刀头/磨头 Shaver Blade Series





目录

CONTENTS


- 刨削刀头 介绍
Shaver Blade Introduction
- 天鹰刨削刀头 特点
Characteristic for Shaver Blade
- 天鹰 优势
Advantage for Aquila
- 刨削刀头 的应用
Application for Shaver System

Aquada Medical

PART 1

刨削刀头 介绍

Shaver Blade Introduction



直径(Dia.):2.0MM-5.8MM, 长度(Length):80MM-133MM

刨刀 介绍

Shaver Blade Introduction

刨刀通常用于软组织切割。在动力系统低速情况下，能产生适宜的负压将组织引入外管窗口中。

天鹰目前刨刀可以适配多种接头，提供长度80MM-300MM, 直径2.0MM-5.3MM选择。

Shaver blades are used for soft tissue cutting. At the low speed of the shaver system, the pressure can be generated to introduce the tissue into the outer tube window.

Aquila shaver blades can be adapted to a variety of joints, providing length 80MM-300MM, diameter 2.0MM-5.3MM.






T=8



T=12



Diamond Blade



Dia.:2.8MM-4.5MM, Length:95MM-355MM

磨头 介绍

Shaver Burr Introduction

磨头通常用于骨质的去除。不同类型的磨头可同时迅速切割软组织和骨组织，提高手术效率；常见的刨削磨头刀刃通常分为6刃、8刃和12刃，刀刃的数量越少，切割性能就越强。

天鹰目前的磨头有8刃、12刃、金刚砂磨头可以选，直径可为2.9MM-5.5MM，长度85.6MM-340MM。

Shaver burr is usually used for bone removal. Different types can quickly cut soft tissue and bone tissue at the same time; Shaver burrs are usually divided into 6, 8, 12 blades, and the fewer the number of blades, the stronger the cutting performance.

PART 2

刨削刀头 特点

Characteristic for Shaver Blade

刨削刀头特点 Characteristic for Shave Blade



AQUILA MEDICAL

A

多种接头适配

Multiple fitting

Dyonics/Arthrex/Linvatec/Stryker/Aquila/Joimax



金属/塑料接头


Metal/plastic connectors



塑料耐高速运转

Plastic resistant to high speed operation

刨削刀头特点 Characteristic for Shave Blade



球头型

Round Burr

鳄齿型

Aggressive
Full Radius

圆柱型

Oval Burr

强力型

Plus

瓜子型

Double
Aggressive Full
Radius

月牙型

Full Radius



AQUILA MEDICAL

B

多种齿形，满足手术不同需求

Multi shapes to meet the different surgery requirement

刨削刀头特点 Characteristic for Shave Blade



B 多种磨头，适用不同场景

A variety of grinding heads, suitable for different scenarios



橄榄型
Oval Burr

球头型
Round Burr

金刚砂
Diamond

万向球型
Multi Angle


磨头耐磨不掉砂
Does not wear off sand

刨削刀头特点 Characteristic for Shave Blade



B 多种角度，更加灵巧、精准定位
A variety of angles, more dexterous,
accurate positioning

万向型磨头 Universal Diamond Burr



耳鼻喉系列 E.N.T (12°~120°)




刨削刀头特点 Characteristic for Shave Blade




C

多种规格型号选择

A variety of specifications to choose



刨刀 Shaver Blade:
直径 (Dia.):2.0MM-5.8MM
长度(Length):84MM-133MM



磨头 Shaver Burr:
直径 (Dia.):2.8MM-4.5MM
长度 (Length):95MM-355MM

刨削刀头特点 Characteristic for Shave Blade



AQUILA MEDICAL

D 适用不同科系

Applicable to different series

- **适配系统:**
Medtronic
- **多角度:**
12°~120°
- **多种规格:**
φ2.8-4.5
- **多种工作长度:**
95-355mm
- **头部选择:**
带刃磨头/金刚砂磨头/刨削刀头



脊柱 Spinal 妇科 Gynaecology 泌尿 Urology

- **适配系统:**
Dyonics/Arthrex/Linvatec/Stryker/Aquila/Joimax
- **多种规格:**
φ3.5-4.4
- **多种工作长度:**
275-340mm
- **刀头选择:**
金刚砂磨头/球头磨头/橄榄型磨头/万向磨头/鳄齿型/月牙型



耳鼻喉
E.N.T.

关节镜(痛风治疗)
Arthroscopy

- **适配系统:**
Dyonics/Arthrex/Linvatec/Stryker/Aquila
- **多种规格:**
φ2.0-5.8
- **多种工作长度:**
84-133mm
- **刀头选择:**
月牙型/瓜子型/强力型/鳄齿型/球头磨头/圆柱型磨头

PART 3

天鹰 优势

Advantage for Aquila

质量优
Good Quality

适配高
High Fitting

规格全
Whole Specification

应用广
More Application

天鹰 优势 Advantage for Aquila

1. 刨刀提供强劲、精细切割，可高效清除软组织；磨头磨削效果佳，刨削刀头转速适用范围:0-12000rpm。
Application range of cutting head speed :0-12000rpm.

2. 内外管直线度严谨，摩擦小，发热低。
The inner and outer tubes have strict straightness, low friction and low heat.

3. 不同规格型号。
Different specifications model.

4. 进口金属材质，特殊加工工艺，不易生锈，多种刨刀/磨头设计，满足不同选择。
Imported metal material

5. 适配多种手柄 Imported metal material
(Dyonics/Arthrex/Linvatec/Stryker/Aquila/Joimax)。

6. 多科系可用（妇科、泌尿科、耳鼻喉、脊柱、关节镜等）。
Multi-disciplinary available

Aquila Medical
天鹰医疗

天鹰 优势

Advantage for Aquila






PART 4

刨削刀头 应用

Application for Shaver Blade

刨削刀头应用

Application for Shaver Blade





AQUILA MEDICAL

天鹰医疗 产品介绍

谢谢观看

Thank You

